


HALIOSEEK® PD-L1/CD8

Data presented at SITC 2017

Halioseek® PD-L1/CD8 is a standardized dual-staining IVD assay which, in addition to PD-L1 detection, provides critical information on the tumor immune infiltrate through the detection of CD8+ cells on the same tissue section. Halioseek® PD-L1/CD8 includes a Digital Pathology (DP) analysis module to determine PD-L1+ and CD8+ cell density as well as the proximity between these cell types.


HALIOSEEK® PD-L1/CD8 DIGITAL PATHOLOGY TOOL PERFORMANCE

CD8+ and PD-L1+ cell detection



- PD-L1+ cell nucleus detected by Halioseek® software
- CD8+ cell detected by Halioseek® software




Proximity index calculation between CD8+ and PD-L1+ cells



- 20µm radius circles surrounding PD-L1+ cells
- PD-L1+ cell with at least 1 CD8+ cell at less than 20µm

Detection Accuracy

34 zones from 25 NSCLC samples (CD8) and 22 zones from 11 NSCLC samples (PD-L1) were analyzed: pathologist absolute cell count of positive cells is highly correlated to DP cell count for both CD8 and PD-L1.




Precision of the proximity index

21 to 28 consecutive slides from 4 NSCLC samples were stained across 4 IHC runs with 3 Halioseek® antibody batches, 3 revelation kit batches on 2 Benchmark XT instruments.

SAMPLES STRATIFICATION WITH DIGITAL PATHOLOGY OUTPUTS


The distribution of PD-L1+ TC % and CD8+ cell densities across 108 NSCLC samples.

PD-L1+ TC % was assessed by pathologists. CD8+ cell density within the tumor area was determined by DP. Different T cell infiltration levels are observed for NSCLC samples with the same PD-L1 status.



Examples of distinct proximity indexes for samples with similar PD-L1 and CD8 levels:

proximity indexes could allow further stratification of patients.



CONCLUSION

Haliousek® PD-L1/CD8 is a new robust IVD assay leveraging the advantages of DP to combine tumor-infiltrating lymphocytes (TILs) and PD-L1 quantification within the tumor microenvironment.

Haliousek® PD-L1/CD8 could have a higher predictive performance than existing IVD tests and could fill a major gap in the management of ICI administration. In a next step we intend to investigate the predictive value of the assay on samples from ICI treated patients.

Haliousek® PD-L1/CD8 is available as a CE-IVD assay and as an RUO* solution in the Rest of The World.

*intended for Research Use Only. Not for use in diagnostic procedures.

REFERENCES

- Conde E et al. Histopathology. 2017
- De Meulenaere A et al. Oncotarget. 2017
- Hiraoka K et al. Br. J. Cancer. 2006
- Parra ER et al. Clin. Cancer Res. 2016
- Schalper KA et al. J. Natl. Cancer Inst. 2015
- Takahashi K et al. Anticancer Res. 2017
- Teng MWL et al. Cancer Res. 2015
- Tumeh PC et al. Nature. 2014
- Velcheti V et al. Lab. Investig. 2014
- Yu H et al. J. Thorac. Oncol. 2016
- Zhuang X et al. Appl. Immunohistochem. Mol. Morphol. 2010

www.haliouseek.com

HALIODX - Luminy Biotech Entreprises

163 Avenue de Luminy - 13288 Marseille Cedex 9 - FRANCE

Phone +33 (0) 4 91 29 30 90 - Fax +33 (0) 4 91 29 30 99 - Email: infos@haliidx.com